

ESCLMC Algorithm: User Guide

This document provides instructions on how to set up and run the ESCLMC prediction algorithm.

1. Prerequisites

- **Hardware:** A machine with a **GPU** (CUDA-enabled) is required to run this program.
- **Software:** Python 3.9 environment. Please install the necessary dependencies using the provided requirements file:

Bash

```
pip install -r requirements.txt
```

2. Setup and Installation

Before running the prediction, you must ensure all files are in the correct location:

1. **Download Resources:** Download the source code, model weights, and configuration files.
2. **Extract Files:**
 - Unzip the **nnUNet weights** directly into the code's root directory.
 - Unzip the **configuration files** and **model checkpoints** into the code's root directory.

Expected Directory Structure:

After extraction, your folder should look like this:

Plaintext

project_root/

```
|—— nnUNet/           # Extracted nnUNet folder
|—— radiomics_config/  # Extracted config folder
|—— min_loss0.pth      # Model weight file
|—— predict.py         # Main script
|—— requirements.txt
|—— ...
```

3. Usage

To run the prediction on a zip file containing DICOM images, open your terminal and run the following command:

Bash

```
python predict.py test_data.zip
```

Replace test_data.zip with the actual path to your input zip file.